





Slide 1









## Motivation: Thermal Management @ mcl



katrin.fladischer@mcl.at

### **Motivation: Thermal Transient Analysis - Simulation** Simulation is used to create the thermal transient. e.g. LED module ∆T(t) Simulated structure function Cooling as function of time Hot **1**P(t) $\log(C_{th\Sigma})$ late

Comparison of simulation and measurements.

Simulated structure function = measured structure function

#### → Validated system

Adjustment of the module's material stack:

Variation of  $\rho, \lambda, c_p$ 

 $R_{th\Sigma}$ 

#### Size Effect – Thermal conductivity



→ Thermal conductivity of 50 nm thin layer of silicon is significantly lower than that of the bulk: from 150 W/mK to 45 W/mK → Temperature dependency follows different law for bulk and thin layers.

LIU *et al.*: MODELING AND DATA FOR THERMAL CONDUCTIVITY OF SOI LAYERS AT HIGH TEMPERATURE IEEE TRANSACTIONS ON ELECTRON DEVICES, VOL. 53, NO. 8, AUGUST 2006

#### Size Effect & Impurities – Thermal conductivity



 $\rightarrow$  The higher the impurity concentration the lower the thermal conductivity.

impurity concentration  $\uparrow$ 

LIU et al.: MODELING AND DATA FOR THERMAL CONDUCTIVITY OF SOI LAYERS AT HIGH TEMPERATURE IEEE TRANSACTIONS ON ELECTRON DEVICES, VOL. 53, NO. 8, AUGUST 2006

Slide 7



#### **Strain – Thermal conductivity**

125 nm Al-films with average grain size of 50 nm (no thickness effect – 3 times the mean free path)

Strain–thermal conductivity coupling:



50% reduction in thermal conductivity at ~0.25% strain (~175 MPa of stress)

mechanical strain **decreases the electron mean free path** →enhanced scattering at the moving grain boundaries

Lee, H. F.; Kumar, S.; Haque, M. A. Role of Mechanical Strain on Thermal Conductivity of Nanoscale Aluminum Films. Acta Mater. 2010, 58, 6619–6627.

#### <u>Time Domain Thermal Reflectance</u> TDTR



Measurement principle

• Optical reflectivity R changes with the temperature.

#### What can be measured?

 Thermal effusivity, diffusivity, conductivity and interfacial thermal resistance of layers

On which scale?

 10 nm to 20 µm thick layers averaged over probe pulse Ø 25 µm

Oven for sample investigation up to 500°C

#### Scanning <u>Th</u>ermal <u>M</u>icroscopy SThM



www.parkafm.com

What can be measured?

- Temperature measurement
- Thermal conductivity
- Topography

On which scale?

Nanometer scale

Slide 9





## Scanning Thermal Microscopy - SThM

#### **SThM**

#### SThM @ MCL: August 2018

#### **Temperatures and thermal conductivities**

#### in the sub-100 nm regime

- lateral (in-plane) resolution
  - < 30 nm



www.parkafm.com



### SPM & SEM combination



#### **SThM: Working Principle**



SEM images of nanolithographed SThM probe from Kelvin Nanotechnology

- **Resistive probe** incorporated in Wheatstone bridge
- Frequenzy-modulated
  measurements
  - $\rightarrow$  3 $\omega$ -method (Fiege et al) Quantitative results

G. Fiege, A. Altes, B. Heiderhoff, L.J. Balk, Quantitative thermal conductivity measurements with nanometre resolution, J. Phys. D-Applied Phys. 32 (1999)

- Cantilever is made of SiN with a thin-film metal wire. Highest resistance of the wire is near the apex of the tip.
- Electrical resistance of thin-film resistor at probe tip correlates with temperature.



General layout of a SThM AFM-based system.

Gomès S., Assy A., and Chapuis P.-O., "Scanning Thermal Microscopy: a review", 2015, Physica Status Solidi (a) 212

Slide 12



#### **SThM: Operating modes**

#### Thermometry

During a scan, heat flows from the **hot sample** to the probe and

changes  $R_p$ 





http://www.powerguru.org/heat-transfer-in-powersemiconductor-devices/

#### **Thermal conductivity**

Hot tip apex acts as nanoscale heat

source





SThM image, 5x5  $\mu m^2$ , of Cu-SiO2-Si stack system, scan in vacuum.

#### **SThM: Thermometry**

#### - Passive mode

- In this mode, a very small electrical current is passed through the probe
- Results in minimal Joule self-heating and enables measurement of electrical resistance R<sub>p</sub> of the probe
- During a scan, heat flows from the hot sample to the probe and changes R<sub>p</sub>



http://www.powerguru.org/heat-transfer-in-power-semiconductor-devices/



#### SThM: Thermal conductivity

- Active mode
- A larger electrical current is passed through the probe, resulting in a significant Joule heating

- Hot tip apex acts as nanoscale heat source
- Thermal conductivity of sample affects SThM probe





SThM image,  $5x5 \ \mu m^2$ , of Cu-SiO2-Si stack system, scan in vacuum.











#### Tungsten (W) – thin film

#### Thermal conductivities of thin films/ interfaces

 Heat managment in metal oxide gas sensors





#### Miniature hotplate: optical image (left) and SEM image (right)

Sideview: Tungsten heating coil for next generation hotplate



SEM image on cross section of CVD grown 500 nm thick W on a 100 nm TiN layer on (100) Si

#### TDTR: Tungsten (W) – thin film – temperature dependent

#### Thermal conductivities ( $\lambda$ ) of thin films







#### x 4 heating cycles

katrin.fladischer@mcl.at

#### TDTR: Tungsten (W) – thin film – temperature dependent

#### Specific heat capacity $(c_p)$ of thin films





#### x 4 heating cycles

0 katrin.fladischer@mcl.at



Thermal conductivities of thin films/ interfaces

#### W\_native = as deposited



W\_long: after a long heat cycle x 4 heating cycles up to 425°C









#### Thermal conductivities of thin films/ interfaces







## W thin film investigated by SThM

#### Thermal conductivity of tungsten thin layers

Sample preparation



SThM on suface



SThM on cross-section

## SThM on cross-section to avoid

- Influence of topography
- Oxidation of surface
- Pollution of surface layer



Ion Slicing to get smooth and even surface



Slide 24 katrin.fladischer@mcl.at

#### Thermal conductivity of tungsten thin layers

Sample preparation



Tungsten sample on sample holder





SThM on cross-section of tungsten sample

#### Thermal conductivity of tungsten thin layers

Topography on cross-section: 500 nm W

100 nm TiN (100) Si substrate



SPM-probe: Arrow-NCR

Slide 26



#### Thermal conductivity of tungsten thin layers

SThM on cross-section

Heating voltage 350 mV, in vacuum:  $1.1*10^{-6}$  mbar

#### W\_native



(a) SThM image

(b) & (c) areas utilized for data evaluation of W and Si at different frequencies

Thermal conductivity, literature values for bulk:

Si ... 149 W/mK W ... 173 W/mK TiN ... 19 W/mK But: TiN layer not visible

$$\lambda_{W_{native}} = 154.2 \pm 4 \text{ W/mK}$$

#### Thermal conductivity of tungsten thin layers

SThM on cross-section

Heating voltage 350 mV, in vacuum:  $1.1*10^{-6}$  mbar

#### W\_short



(a) SThM image

(b) & (c) areas utilized for data evaluation of W and Si at different frequencies

$$\lambda_{W\_short} = 151.8 \pm 4 \text{ W/mK}$$

#### Thermal conductivity of tungsten thin layers

SThM on cross-section

Heating voltage 350 mV, in vacuum:  $1.1*10^{-6}$  mbar

#### W\_long



(a) SThM image

(b) & (c) areas utilized for data evaluation of W and Si at different frequencies

$$\lambda_{W_{long}} = 155.6 \pm 4 \text{ W/mK}$$

#### Tungsten (W) – thin film – at 30°C SPM topview – Grain sizes



#### Tungsten (W) – thin film – at 30°C

#### **SPM** topview – Grain sizes



D. Choi, The electron scattering at grain boundaries in tungsten films, Microelectron. Eng. 122 (2014) 5–8. doi:10.1016/j.mee.2014.03.012:

Grain boundary scattering  $\rightarrow$  increase of electrical resistivity due to

- surface scattering (Fuchs-Sonderheimer model  $\Delta \sigma_{FS}$ ) and
- grain boundary scattering (Mayadas-Shatzkes model  $\Delta \sigma_{MS}$ )
- → Caclulation of  $\lambda_{calc}$  out of grain size

Sample	Average grain	$\lambda_{calc}$	$\lambda_{TDTR}$	λ <sub>sthm</sub>
	size [nm]	[W/mK]	[W/mK]	[W/mK]
W_native	136 ± 51	$152^{+5}_{-11}$	157 ± 8	153.6 ±4
W_short	105 ± 50	$146^{+8}_{-19}$	149.8 ± 0.5	151.8 ±4
W_long	177 ± 84	$156^{+5}_{-13}$	195 ± 18	155.6 ±4

Slide 31



GIS

GEN

ople Holder

SPM

Sample



## Outlook

Outlook

#### Outlook – SThM at MCL







Simulation cell representing layered system. Gibbons, Bebek, Kang, Stanley, and Estreicher, J. of Applied Physics 118 (2015).

Topography showing steps of strontium titanate: Image size 1.1µm https://www.nanosurf.com/en/application/547-topography-strontium-titanate-flexafm

#### SThM - SEM characterization

- Study thermal transport mechanism, phonons /electrons?
- In situ loading of hydrogen etc.



Materials Center Leoben Forschung GmbH www.mcl.at

# Any questions?

## Thank you for your attention!

#### Acknowledgement:

COMÉT

Competence Centers for Excellent Technologies bm

This project has received funding from the Electronic Component Systems for European Leadership Joint Undertaking under grant agreement No 692480 (IoSense). This Joint Undertaking receives support from the European Union's Horizon 2020 research and innovation programme and Germany, Spain, Austria, Belgium, Slovakia.

Das Land



